

AL04-06-230

Ultrasonic Water Bath 6 Liter (heating)

*“Reliable easy to use equipment for finest cleaning applications.
Availability of cover constructed basket holder keeping a clean bench.
Must have cleaning solution in any lab.”*



APPLICATIONS:

Suitable for application listed for AL04-04-230 but with larger basket dimensions and deeper tank for long instruments or a larger amount of tools or cleaning items for:

- The use in the laboratory for cleaning and degassing
- The processing of drills and burs in the **dental practice**
- Cleaning of small items in the **medical sector**
- Cleaning of small mechanic item and tools in the **industry**

FEATURES:

- **Transducer tank** made of **cavitation-proof stainless steel**.
- Housing made of stainless steel, easy to keep hygienically clean.
- Sandwich performance transducer systems.
- **Degas function** for the efficient degassing of the cleaning liquid and for special laboratory applications.
- Auto-Degas function for an automatic degassing cycle, e.g. after freshly mixed cleaning liquid has been filled.
- Indication of remaining cleaning time via diodes.
- Splashwater-proof operating panel.
- **Automatic switch-off** after 12 h operation (if no buttons have been activated in this period) to prevent unintended permanent operation.
- Electronic turning knobs.

TECHNICAL SPECIFICATION

Size (Liter)	6l
Tank max. volume (app. Liter)	5.75
Tank working volume (app. Liter)	4.3
Tank inside dimensions	300 x 151 x 150 (B x T x H / ca. mm)
Unite outside dimensions	365 x 186 x 264 (B x T x H / ca. mm)
Ultrasonic peak max. (W)	600 Watt
Ultrasonic power effective (W)	150 Watt
Total power consumption (W)	550 Watt
Heating Power	400 Watt
Frequency (kHz)	38 kHz
Power supplier (Vac)	100-120 / 220-240 (Vac)
Weight (app. Kg)	5.1
Material outside (stainless steel)	V2A 1.4301 (AISI 304)
Material tank (stainless steel)	V2A 1.4301 (AISI 304)
Tank thickness (ca. mm)	0.8

ACCESSORIES

- Basket inside dimensions (ca. mm): 255 x 115 x 75
- Basket volume (app. Kg): 5.0

Code	Description
ALo4-06-151	Basket
ALo4-06-150	Plastic cover
ALo4-06-152	Set 2 breakers incl. covers
ALo4-01-153	Cover for beaker
ALo4-01-154	Rubber ring for beaker
ALo4-01-155	Immersion basket (stainless steel)
ALo4-01-156	Jewelry Rack



CHEMICALS

Code	Description
D-cleaner range	Cleaning chemicals for dental applications
T-cleaner range	Cleaning chemicals for industrial applications
L-cleaner range	Cleaning chemicals for laboratory applications
O-cleaner range	Cleaning chemicals for the optician

*Prices and cleaning information can be found on www.advantage-lab.com