

Accessories Vertical Freezers

Upright Drawer Rack

Strong stainless steel racks produced with drawers for easy practical access to your boxes of samples



Upright Side Rack

Stainless steel rack with individual compartments for each box. Very solid and strong. On this model the box is taken out from the side.



Model	Description	Box size (mm)	No. of box – Depth	No. of box – Height	Total	Dimension WxDxH (mm)
AL07-99-220	Upright Drawer Rack 50-34	50	3	4	12	141 x 434 x 231
AL07-99-221	Upright Drawer Rack 50-35	50	3	5	15	141 x 434 x 290
AL07-99-222	Upright Drawer Rack 50-44	50	4	4	16	141 x 567 x 231
AL07-99-223	Upright Drawer Rack 50-45	50	4	5	20	141 x 567 x 290
AL07-99-224	Upright Drawer Rack 75-33	75	3	3	9	141 x 434 x 244
AL07-99-225	Upright Drawer Rack 75-43	75	4	3	12	141 x 567 x 244
AL07-99-226	Upright Drawer Rack 96-32	96	3	2	6	141 x 434 x 204
AL07-99-227	Upright Drawer Rack 96-33	96	3	3	9	141 x 434 x 306
AL07-99-228	Upright Drawer Rack 96-42	96	4	2	8	141 x 567 x 204
AL07-99-229	Upright Drawer Rack 96-43	96	4	3	12	141 x 567 x 306
AL07-99-250	Upright Side Rack 50-34	50	3	4	12	139 x 414 x 224
AL07-99-251	Upright Side Rack 50-35	50	3	5	15	139 x 414 x 279
AL07-99-252	Upright Side Rack 50-44	50	4	4	16	139 x 549 x 224
AL07-99-253	Upright Side Rack 50-45	50	4	5	20	139 x 549 x 279
AL07-99-254	Upright Side Rack 75-33	75	3	3	9	139 x 414 x 240
AL07-99-255	Upright Side Rack 75-43	75	4	3	12	139 x 549 x 240
AL07-99-256	Upright Side Rack 96-32	96	3	2	6	139 x 414 x 202
AL07-99-257	Upright Side Rack 96-33	96	3	3	9	139 x 414 x 301
AL07-99-258	Upright Side Rack 96-42	96	4	2	8	139 x 549 x 202
AL07-99-259	Upright Side Rack 96-43	96	4	3	12	139 x 549 x 301

Freezer Boxes

There are many different plastic boxes on the market. Our chosen manufacturer is a European based company with many years of experience in making strong and easy to use storage boxes of a strong polypropylene mixture. The box is equipped with a standard 9x9 divider to provide a storage capacity of 81 tubes with a diameter up to 12,5 mm in each box.



Model	Box size (mm)	Divider	Total no. of vials	Material	Dimensions (mm)
AL07-99-150	50	9 x 9	81	Polypropylene	130 x 130 x 50
AL07-99-151	75	9 x 9	81	Polypropylene	130 x 130 x 75
AL07-99-152	96	9 x 9	81	Polypropylene	130 x 130 x 96

AL07-23-230 Freezer Upright, ultra low -40 / -86°C, 826L

Options	Description	Box size (mm)	No. of boxes per Rack - Dept x Height	Total no. of boxes	Dimension WxDxH (mm)
1	24 x Upright drawer Rack 50-44 +6 x Upright drawer Rack 50-45 (=30 racks)	50	4x4 4x5	384 + 120 (=504 boxes)	141 x 567 x 231 141 x 567 x 290
2	30 x Upright drawer Rack 75-43	75	4x3	360	141 x 567 x 244
3	24 x Upright drawer Rack 96-42 +6 x Upright drawer Rack 96-43 (=30 racks)	96	4x2 4x3	192 + 72 (=264 boxes)	141 x 567 x 204 141 x 567 x 306
4	24 x Upright Side Rack 50-44 +6 x Upright Side Rack 50-45 (= 30 racks)	50	4x4 4x5	384 + 120 (=504 boxes)	139 x 549 x 224 139 x 549 x 279
5	24 x Upright Side Rack 75-43 +6 x Upright Side Rack 75-44 (= 29 racks)	75	4x3 4x4	288 + 80 (=368 boxes)	139 x 549 x 240 139 x 549 x 324
6	24 x Upright Side Rack 96-42 +5 x Upright Side Rack 96-43 (= 29 racks)	96	4x2 4x3	192 + 60 (=252 boxes)	139 x 549 x 202 139 x 549 x 301

AL07-21-230 Freezer Upright, ultra low -40 / -86°C, 680L

Options	Description	Box size (mm)	No. of boxes per Rack - Dept x Height	Total no. of boxes	Dimension WxDxH (mm)
1	20 x Upright Drawer Rack 50-44 +5 x Upright Drawer Rack 50-45 (=25 racks)	50	4x4 4x5	320 + 100 (=420 boxes)	141 x 567 x 231 141 x 567 x 290
2	25 x Upright drawer Rack 75-43	75	4x3	300	141 x 567 x 244
3	20 x Upright Drawer Rack 96-42 +5 x Upright Drawer Rack 96-43 (=25 racks)	96	4x2 4x3	160 + 60 (=220 boxes)	141 x 567 x 204 141 x 567 x 306
4	20 x Upright Side Rack 50-44 +5 x Upright Side Rack 50-45 (= 25 racks)	50	4x4 4x5	320 + 100 (=420 boxes)	139 x 549 x 224 139 x 549 x 279
5	25 x Upright Side Rack 75-43	75	4x3	300	139 x 549 x 240
6	20 x Upright Side Rack 96-42 +5 x Upright Side Rack 96-43 (= 25 racks)	96	4x2 4x3	160 + 60 (=220 boxes)	139 x 549 x 202 139 x 549 x 301

AL07-19-230 Freezer Upright, ultra low -40 / -86°C, 585L

Options	Description	Box size (mm)	No. of boxes per Rack - Dept x Height	Total no. of boxes	Dimension WxDxH (mm)
1	20 x Upright Drawer Rack 50-44	50	4x4	320	141 x 567 x 231
2	20 x Upright Drawer Rack 75-43	75	4x3	240	141 x 567 x 244
3	20 x Upright Drawer Rack 96-42	96	4x2	160	141 x 567 x 204
4	20 x Upright Side Rack 50-45	50	4x5	400	139 x 549 x 279
5	20 x Upright Side Rack 75-43	75	4x3	240	139 x 549 x 240
6	20 x Upright Side Rack 96-42	96	4x2	160	139 x 549 x 202

AL07-13-230 Freezer Upright, ultra low -40 / -86°C, 413L

AL07-12-230 Freezer Upright, low -20 / -40°C, 413L

Options	Description	Box size (mm)	No. of boxes per Rack - Dept x Height	Total no. of boxes	Dimension WxDxH (mm)
1	12 x Upright Drawer Rack 50-44 +3 x Upright Drawer Rack 50-45 (=15 racks)	50	4 x 4 4 x 5	192 +60 (=252 boxes)	141 x 567 x 231 141 x 567 x 290
2	15 x Upright Drawer Rack 75-43	75	4 x 3	180	141 x 567 x 244
3	15 x Upright Drawer Rack 96-42	96	4 x 2	120	141 x 567 x 204
4	12 x Upright Side Rack 50-44 +3 x Upright Side Rack 50-45 (=15 racks)	50	4 x 4 4 x 5	192 +60 (=252 boxes)	139 x 549 x 224 139 x 549 x 279
5	15 x Upright Side Rack 75-43	75	4 x 3	180	139 x 549 x 240
6	12 x Upright Side Rack 50-44 +3 x Upright Side Rack 50-45 (=15 racks)	96	4 x 2 4 x 3	96 +36 (=132 boxes)	139 x 549 x 202 139 x 549 x 301

AL07-15-230 Freezer Upright, ultra low -40 / -86°C, 115L

AL07-14-230 Freezer Upright, low -20 / -40°C, 115L

Options	Description	Box size (mm)	No. of boxes per Rack - Dept x Height	Total no. of boxes	Dimension WxDxH (mm)
1	4 x Upright Drawer Rack 50-34	50	3 x 4	48	141 x 434 x 231
2	4 x Upright Drawer Rack 75-33	75	3 x 3	136	141 x 434 x 244
3	4 x Upright Drawer Rack 96-32	96	3 x 2	24	141 x 434 x 204
4	4 x Upright Side Rack 50-34	50	3 x 4	48	139 x 414 x 224
5	4 x Upright Side Rack 75-33	75	3 x 3	36	139 x 414 x 240
6	4 x Upright Side Rack 96-32	96	3 x 2	24	139 x 414 x 202

AL07-17-230 Freezer Upright, ultra low -40 / -86°C, 55L

AL07-16-230 Freezer Upright, low -20 / -40°C, 55L

Options	Description	Box size (mm)	No. of boxes per Rack - Dept x Height	Total no. of boxes	Dimension WxDxH (mm)
1	2 x Upright Drawer Rack 50-35	50	3 x 5	30	141 x 434 x 290
2	2 x Upright Drawer Rack 75-33	75	3 x 3	18	141 x 434 x 244
3	2 x Upright Drawer Rack 96-33	96	3 x 3	18	141 x 434 x 306
4	2 x Upright Side Rack 50-35	50	3 x 5	30	139 x 414 x 279
5	2 x Upright Side Rack 75-33	75	3 x 3	18	139 x 414 x 240
6	2 x Upright Side Rack 96-33	96	3 x 3	18	139 x 414 x 301